

11

10

9

8

7

6

5

4

3

2

1

DASH NO

TYPE OR MODEL

NEXT ASSEMBLY

-

SAPPHIRE

TBD

REF DES

QTY

PART NUMBER

DESCRIPTION OR MATERIAL

ITEM

REF

P1060967

SCH, STAND CONTROLLER

1

2

3

4

1

P1060969003

DBB, STAND CONTROLLER

5

1

P1060971004

BOM, STAND CONTROLLER

6

7

8

9

10

4

KEYING TABLE

REF DES

ADD KEYING PLUG SHOWN IN BOM TO CONTACT NO.

J22

5

J9

7

J29

12

J30

12

J28

13

J7

6

J3

7

J4

7

J12

17

J6

1

REVISION BLOCK

DESC.

DWG OR PART NO.

REV

SCHEM

P1060967

D

DETAIL

P1060969003

C

BOM

P1060971004

D

NOTES: UNLESS OTHERWISE SPECIFIED.

1. ASSEMBLY WORKMANSHIP TO BE IN ACCORDANCE TO IPC-A-610 CLASS 2, LATEST REVISION.

2. FINISHED ASSEMBLY MUST BE RoHS COMPLIANT.

3 MARK DASH NUMBER AND REVISION OF ASSEMBLY, REVISION ONLY OF SCHEMATIC, AND COUNTRY OF ASSEMBLY. IN LOCATIONS SHOWN. MARKING SHALL BE PERMANENT AND REMAIN INTACT AND LEGIBLE THROUGHOUT ALL PROCESSES.

4 INSTALL KEYING PLUGS REFERENCED IN BILL OF MATERIALS PER KEYING TABLE.

THIS DOCUMENT CONTAINS INFORMATION WHICH IS THE PROPRIETARY PROPERTY OF VARIAN MEDICAL SYSTEMS. REPRODUCTION, DISCLOSURE, OR RELEASE TO OTHERS, MANUALLY OR ELECTRONICALLY, WITHOUT THE PRIOR WRITTEN CONSENT OF VARIAN MEDICAL SYSTEMS IS STRICTLY PROHIBITED.

DO NOT SCALE DRAWING

DESCRIPTION OF CHANGE

DISC. 003. ADDED 004
CHANGED SCH TO REV D
CHANGED DETAIL TO 003 REV C
CHANGED BOM TO 004 REV D
REMOVED DIGITAL DATA TABLE
CHANGED TOP AND BOTTOM VIEWS
REMOVED NOTES 5 AND 6
REMOVED REF TO REWORK DOC

APVD

ECO

DFT

CHK

DATE

REV

-

-

KG

-

2024-10-31

D

UNLESS OTHERWISE SPEC:
DIMENSIONS ARE IN INCHES
ANGS=

FIN ✓

D

P1060966

D

SIZE

DRAWING NUMBER

REV

ACAD=60.

PCB ASSY,
STAND CONTROLLER

P1060966004

PCB ASSY, STAND CONTROLLER

P1060966003

DISCONTINUED

P1060966002

DISCONTINUED

P1060966001

DISCONTINUED

DASH

DESCRIPTION

varian

DRAWN

DATE

APPROVED

DATE

SCALE

K.GREENBERG

2023 MAY 24

-

-

2/1

CHECKED

DATE

APPROVED

DATE

SHEET

-

-

-

1 OF 3

11

10

9

8

7

6

5

4

3

2

1

11

10

9

8

7

6

5

4

3

2

1


CM# 000060021312

Released

Ver: 03

Rev: D

11. November 2024



STAND CONTROLLER
DETAIL: P1060969003

ASSY: P1060966
SCH: P1060967

TOP SIDE VIEW

varian				SCALE 2/1	
				SHEET 2 OF 3	
PCB ASSY, STAND CONTROLLER					
UNLESS OTHERWISE SPEC:					
ANG±=		-	FIN		✓
X±=		.1	XX±=		.03
			XXX±=		.005
D		P1060966		D	
SIZE		DRAWING NUMBER			REV

CM# 000060021312

Released
Ver: 03

Rev: D

11. November 2024

Signature File



Document: P1060966/ONC/000/03
Description: PCB ASSY, STAND CONTROLLER

Change Master Number: 60021312
ECN Number: 200109931
ECN Description:

Signature list

Date	Time	User	Status
11, November 2024	18:43:59	Gary Yen	Verification Rel
06, November 2024	01:35:39	Ge Song	Approved
06, November 2024	00:29:29	Kevin Greenberg	In Eng Approval
06, November 2024	00:03:41	Gary Yen	In Checking
05, November 2024	18:44:53	Gary Yen	In Drafting
04, November 2024	19:26:59	Gary Yen	Unreleased/Orig

Signature File



Document: P1060966/ONC/000/03
Description: PCB ASSY, STAND CONTROLLER

Change Master Number: 60021312
ECN Number: 200109931
ECN Description:

Signature list

Date	Time	User	Status
11, November 2024	18:45:05	Gary Yen	Released
11, November 2024	18:43:59	Gary Yen	Verification Rel
06, November 2024	01:35:39	Ge Song	Approved
06, November 2024	00:29:29	Kevin Greenberg	In Eng Approval
06, November 2024	00:03:41	Gary Yen	In Checking
05, November 2024	18:44:53	Gary Yen	In Drafting
04, November 2024	19:26:59	Gary Yen	Unreleased/Orig